







Pride Mobility Products is a world leader in the design, development and manufacture of mobility products - powerchairs, scooters and lift chairs - for people with disabilities and mobility impairments. With flagship product ranges such as the Jazzy Range, the British manufactured Hudson Range and our innovative Zero Turn technology – Pride Mobility are committed to helping improving people's lives through quality, innovative products.





Red

White

Blue

Yellow

Quartz

Orange

One-hand feather-touch disassembly



21stone 61bs 3.7mph

Weight Capacity

Maximum Speed



Turning Radius

WEIGHT CAPACITY

21 stone 6 lbs (136 kg)

MAX SPEED

3.7 mph (5.9 kmh)

TURNING RADIUS

645 mm (25.5")

OVERALL LENGTH

840 mm (33")

OVERALL WIDTH

574 mm (22.5")

DRIVE WHEELS

229 mm (solid tyres)

BATTERIES REQUIRED

(2) 18 AH (deep cycle)



18 stone 8 lbs (120 kg)

MAX SPEED

3.7 mph (5.9 kmh)

TURNING RADIUS

900 mm (35.4")

OVERALL LENGTH

880 mm (34.6")

OVERALL WIDTH

600 mm (23.6")

DRIVE WHEELS

216 mm (solid tyres)

BATTERIES REQUIRED

(2) 6.6 AH (Li-ion battery)









17 stone 12 lbs 3.6 mph

Weight Capacity

Turning Radius

WEIGHT CAPACITY

17 stone 12 lbs (113 kg)

MAX SPEED

3.6 mph (5.8 kmh)

TURNING RADIUS

616 mm (24.25")

OVERALL LENGTH

1016 mm (40")

OVERALL WIDTH

597 mm (23.5")

DRIVE WHEELS

305 mm (pneumatic)

BATTERIES REQUIRED

(1) 18.2 AH (Li-ion battery)



18 stone 8 lbs (120 kg)

MAX SPEED

4 mph (6.4 kmh)

TURNING RADIUS

650 mm (25.59")

OVERALL LENGTH

1100 mm (43.3")

OVERALL WIDTH

620 mm (24.4")

DRIVE WHEELS

250 mm (solid tyres)

BATTERIES REQUIRED

(1) 50 AH (Li-ion battery)





Turning Radius



Maximum Speed

21stone 61bs

Weight Capacity



WEIGHT CAPACITY

21 stone 6 lbs (136 kg)

MAX SPEED

4 mph (6.4 kmh)

TURNING RADIUS

551 mm (22.5")

OVERALL LENGTH

1118 mm (44")

OVERALL WIDTH

603 mm (23.75")

DRIVE WHEELS

254 mm (solid tyres)

BATTERIES REQUIRED

(2) 30 AH (deep cycle)



21 stone 6 lbs (136 kg)

MAX SPEED

4 mph (6.4 kmh)

TURNING RADIUS

546 mm (21.5")

OVERALL LENGTH

1099 mm (43.25")

OVERALL WIDTH

603 mm (23.75")

DRIVE WHEELS

254 mm (solid tyres)

BATTERIES REQUIRED

(2) 30 AH (deep cycle)





Turning Radius



Up To

Hmph

Maximum Speed

21stone 61bs

Weight Capacity



WEIGHT CAPACITY

21 stone 6 lbs (136 kg)

MAX SPEED

4 mph (6.4 kmh)

TURNING RADIUS

521 mm (20.5")

OVERALL LENGTH

1022 mm (40.25")

OVERALL WIDTH

622 mm (24.5")

DRIVE WHEELS

355 mm (solid tyres)

BATTERIES REQUIRED

(2) 50 AH (deep cycle)

Three optional seat sizes:

- 1. 16"W Adjustable depth of 16"-18"
- 2. 18"W Adjustable depth of 18"-20"
- 3. 20"W Adjustable depth of 18"-20"



21 stone 6 lbs (136 kg)

MAX SPEED

4 mph (6.4 kmh)

TURNING RADIUS

528 mm (20.75")

OVERALL LENGTH

1067 mm (42")

OVERALL WIDTH

590 mm (23.2")

DRIVE WHEELS

254 mm (solid tyres)

BATTERIES REQUIRED

(2) 40 AH (deep cycle)

Available in three sizes:

1. 41 x 41-46cm

2. 46 x 46-51cm

3. 51 x 51-56cm









3 different

colour panels available

Weight Capacity

Hmph

Maximum Speed

530mm
Turning Radius



WEIGHT CAPACITY

21 stone 6 lbs (136 kg)

MAX SPEED

4 mph (6.4 kmh)

TURNING RADIUS

530 mm (21")

OVERALL LENGTH

890 mm (35")

OVERALL WIDTH

635 mm (25")

DRIVE WHEELS

355 mm (pneumatic)

BATTERIES REQUIRED

(2) 50 AH (deep cycle)

Two configurations available:

- 1. Power tilt only, manual recline and standard swing-away leg rests.
- 2. Power tilt, power recline and manual elevating leg rests.

ACCESSORIES



XLR USB Charger & Phone Mounts

Charge your phone on the go! This can be used while your powerchair or scooter is in motion. Compatible with most mobile phones.



Crutch & Walker Holders

Quick and easy access for your cane, crutch or walking stick with our selection of holders.



Baskets & Bags

Our handy and versatile range of additional large baskets and storage bags make transporting your essentials hassle free.



Foldable Travel Bags

Light and durable travel bags for the Quest, iGo, iGo Plus and iGo Fold.



Rear Oxygen Tank Holder

Available with most models, featuring 2 thumb screws for safe transportation.



Cup Holders

Satisfy your thirst on the move with our handy drinks holders, compatible with most of our range of products.







Also from Pride:



Travel Scooter Range



Large Scooter Range



Lift Chair Range

Pride Mobility LTD

32 Wedgwood Road Tel: 01869 324600

Bicester, Oxon Email: sales@pride-mobility.co.uk

OX26 4UL, England Web: www.pride-mobility.co.uk











Pride Mobility LTD reserve the right to alter, at any time, without prior notice, the design, specification, packaging or price without the incurring any obligation.